

### **Certificate of Analysis**

COA No: CA BMM-0018

Version: 09

## ImmoMix™ Red

For research or further manufacturing use only

Catalog No:	BIO-25022
Lot No:	PM356-B112330
Storage Conditions:	-20°C
Component Lot No:	IMXR-222311A
Expiry date:	December 2024

### **Quality Control Parameters**

Analysis	Specification	Result
Functional	A 0.9Kb fragment is amplified with a dilution series of human genomic DNA, using standard conditions and 30 cycles. Single distinct bands were observed with agarose gel electrophoresis (ethidium stained).	Passed
DNA contamination	Quantitative PCR analysis with no template. Presence of <i>E. coli</i> and mouse genomic DNA checked. Test sample must amplify in line with a reference sample.	Passed
DNase contamination	Incubation of a 1Kb double stranded DNA fragment. Incubation for 4 hours at $37^{\circ}$ C with dilution series of DNase I. Analysed by agarose gel electrophoresis. Test sample must show less degradation than the limit of detection $2.5 \times 10^{-3}$ U DNase.	Passed

QA / QC Representative:



Alberta Newton

Date: 18th November 2022



### **Certificate of Analysis**

COA No: CA\_XBB-0014

Version: 08

# MgCl<sub>2</sub> Solution, 50mM

For research or further manufacturing use only

Catalog No:	BIO-25022
Lot No:	PM356-B112330
Storage Conditions:	-20°C
Component Lot No:	MG-2031.015
Expiry date:	December 2024

### **Quality Control Parameters**

Analysis	Specification	Result
Functional	Fragments of sizes 800bp and 3000bp are amplified with a dilution series of BIOTAQ <sup>TM</sup> DNA Polymerase, using standard conditions and 30 cycles. Single distinct bands were observed with agarose gel electrophoresis (ethidium stained).	Passed
DNA contamination	Quantitative PCR analysis with no template. Presence of <i>E. coli</i> and mouse genomic DNA checked. Test sample must amplify in line with a reference sample.	Passed
DNase contamination	Incubation of a 1Kb double stranded DNA fragment. Incubation for 4 hours at 37°C with dilution series of DNase I. Analysed by agarose gel electrophoresis. Test sample must show less degradation than the limit of detection 2.5 x 10 <sup>-3</sup> U DNase.	Passed

QA / QC Representative:



Alberta Newton

Date: 18<sup>th</sup> November 2022

United Kingdom

Tel: +44 (0)20 8830 5300 Fax: +44 (0)20 8452 2822 <u>USA</u>

Tel: +1 901.382.8716 Fax: +1 901.382.0027 <u>Germany</u>

Tel: +49 (0)3371 60222 00 Fax: +49 (0)3371 60222 01 <u>Australia</u>

Tel: +61 (0)2 9209 4180 Fax: +61 (0)2 9209 4763