

# Reverse Transcriptase

For research or further manufacturing use only

|                     |                |
|---------------------|----------------|
| Catalog No:         | BIO-65043      |
| Lot No:             | RA384-B110300  |
| Storage Conditions: | -20°C          |
| Component Lot No:   | TRT-222108A    |
| Expiry date:        | September 2024 |

## Quality Control Parameters

| Analysis                      | Specification  | Result |
|-------------------------------|--|--------|
| Functional                    | Fragments of sizes 1.2Kb and 6.5Kb were reverse transcribed, using standard conditions. Single distinct bands were observed with agarose gel electrophoresis (ethidium stained). | Passed |
| Endonuclease contamination    | Super coiled DNA plasmid was incubated with the reverse transcriptase for 1 hour at 37°C, the absence of nicking and cutting is shown by agarose gel electrophoresis.            | Passed |
| DNase and RNase contamination | A DNA and RNA fragment were Incubated with the reverse transcriptase for 1 hour at 37°C. < 1% degradation was observed.  | Passed |

QA / QC Representative:



Andrew Galeeba-M

Date: 5<sup>th</sup> September 2022

**United Kingdom**

Tel: +44 (0)20 8830 5300  
Fax: +44 (0)20 8452 2822

**USA**

Tel: +1 901.382.8716  
Fax: +1 901.382.0027

**Germany**

Tel: +49 (0)3371 60222 00  
Fax: +49 (0)3371 60222 01

**Australia**

Tel: +61 (0)2 9209 4180  
Fax: +61 (0)2 9209 4763

**RT Buffer**

For research or further manufacturing use only

|                     |                |
|---------------------|----------------|
| Catalog No:         | BIO-65043      |
| Lot No:             | RA384-B110300  |
| Storage Conditions: | -20°C          |
| Component Lot No:   | TRTB-222108A   |
| Expiry date:        | September 2024 |

**Quality Control Parameters**

| Analysis            | Specification   | Result |
|---------------------|---|--------|
| Functional          | Fragment of size 1Kb was reverse transcribed with BioScript™, with a template dilution series, using standard conditions. Single distinct bands were observed with agarose gel electrophoresis (ethidium stained).  | Passed |
| DNA contamination   | Quantitative PCR analysis with no template. Presence of <i>E. coli</i> and mouse genomic DNA checked. Test sample must amplify in line with a reference sample.   | Passed |
| DNase contamination | Incubation of a 1Kb double stranded DNA fragment. Incubation for 4 hours at 37°C with dilution series of DNase I. Analysed by agarose gel electrophoresis. Test sample must show less degradation than the limit of detection $2.5 \times 10^{-3}$ U DNase. | Passed |
| RNase contamination | Quantitative PCR analysis with high and low RNase standards. Test sample must show less RNase activity than the limit of detection $9.7 \times 10^{-3}$ ng/μL RNase.  | Passed |

QA / QC Representative:



Andrew Galeeba-M

Date: 5<sup>th</sup> September 2022

**United Kingdom**

Tel: +44 (0)20 8830 5300  
Fax: +44 (0)20 8452 2822

**USA**

Tel: +1 901.382.8716  
Fax: +1 901.382.0027

**Germany**

Tel: +49 (0)3371 60222 00  
Fax: +49 (0)3371 60222 01

**Australia**

Tel: +61 (0)2 9209 4180  
Fax: +61 (0)2 9209 4763

**RNase Inhibitor**

Suitable for Research and further Manufacturing Use

|                     |                |
|---------------------|----------------|
| Catalog No:         | BIO-65043      |
| Lot No:             | RA384-B110300  |
| Storage Conditions: | -20°C          |
| Component Lot No:   | RI-122308A     |
| Expiry date:        | September 2024 |

**Quality Control Parameters**

| Analysis   | Specification   | Result |
|------------|---|--------|
| Inhibition | Test level of inhibition by incubating total RNA with concentration gradient of RNase A. Bands were observed with agarose gel electrophoresis (ethidium stained). | Passed |

QA / QC Representative:



Andrew Galeeba-M

Date: 5<sup>th</sup> September 2022

**United Kingdom**

Tel: +44 (0)20 8830 5300  
Fax: +44 (0)20 8452 2822

**USA**

Tel: +1 901.382.8716  
Fax: +1 901.382.0027

**Germany**

Tel: +49 (0)3371 60222 00  
Fax: +49 (0)3371 60222 01

**Australia**

Tel: +61 (0)2 9209 4180  
Fax: +61 (0)2 9209 4763

|  |                     |                |
|--|---------------------|----------------|
| <p><b>dNTP Mix 10mM</b></p> <p>Suitable for Research and further Manufacturing Use</p> | Catalog No:         | BIO-65043      |
|  | Lot No:             | RA384-B110300  |
|  | Storage Conditions: | -20°C          |
|  | Component Lot No:   | DM1-122208A    |
|  | Expiry date:        | September 2024 |

**Quality Control Parameters**

| Analysis            | Specification   | Result |
|---------------------|---|--------|
| Functional          | A 800bp human genomic DNA fragment is amplified with a dilution series of enzymes, using standard conditions and 30 cycles. Single distinct bands were observed with agarose gel electrophoresis (ethidium stained).  | Passed |
| DNA contamination   | Quantitative PCR analysis with no template. Presence of <i>E. coli</i> and mouse genomic DNA checked. Test sample must amplify in line with a reference sample.   | Passed |
| DNase contamination | Incubation of a 1Kb double stranded DNA fragment. Incubation for 4 hours at 37°C with dilution series of DNase I. Analysed by agarose gel electrophoresis. Test sample must show less degradation than the limit of detection $2.5 \times 10^{-3}$ U DNase. | Passed |
| RNase contamination | Quantitative PCR analysis with high and low RNase standards. Test sample must show less RNase activity than the limit of detection $9.7 \times 10^{-3}$ ng/ $\mu$ L RNase.  | Passed |
| Nicking Activity    | Incubation of dNTP Mix with supercoiled control plasmid. Analysed by agarose gel electrophoresis. Test sample does not show an increase of linearized or relaxed plasmid.   | Passed |

Authorised by QA / QC Representative:



Andrew Galeeba-M

Date: 5<sup>th</sup> September 2022

**United Kingdom**

Tel: +44 (0)20 8830 5300  
Fax: +44 (0)20 8452 2822

**USA**

Tel: +1 901.382.8716  
Fax: +1 901.382.0027

**Germany**

Tel: +49 (0)3371 60222 00  
Fax: +49 (0)3371 60222 01

**Australia**

Tel: +61 (0)2 9209 4180  
Fax: +61 (0)2 9209 4763

**Oligo (dT)<sub>18</sub>**

For research or further manufacturing use only

|                     |                |
|---------------------|----------------|
| Catalog No:         | BIO-65043      |
| Lot No:             | RA384-B110300  |
| Storage Conditions: | -20°C          |
| Component Lot No:   | ODT-122108A    |
| Expiry date:        | September 2024 |

**Quality Control Parameters**

| Analysis            | Specification   | Result |
|---------------------|---|--------|
| Functional          | A 1Kb fragment is reverse transcribed from Poly A RNA with a dilution series of BioScript™, using standard conditions. Single distinct bands were observed with agarose gel electrophoresis (ethidium stained). | Passed |
| RNase contamination | Quantitative PCR analysis with high and low RNase standards. Test sample must show less RNase activity than the limit of detection $9.7 \times 10^{-3}$ ng/μL RNase.  | Passed |

QA / QC Representative:



Andrew Galeeba-M

Date: 5<sup>th</sup> September 2022

**United Kingdom**

Tel: +44 (0)20 8830 5300  
Fax: +44 (0)20 8452 2822

**USA**

Tel: +1 901.382.8716  
Fax: +1 901.382.0027

**Germany**

Tel: +49 (0)3371 60222 00  
Fax: +49 (0)3371 60222 01

**Australia**

Tel: +61 (0)2 9209 4180  
Fax: +61 (0)2 9209 4763

## Random Hexamer Primer

For research or further manufacturing use only

|                     |                |
|---------------------|----------------|
| Catalog No:         | BIO-65043      |
| Lot No:             | RA384-B110300  |
| Storage Conditions: | -20°C          |
| Component Lot No:   | RHP-122108A    |
| Expiry date:        | September 2024 |

### Quality Control Parameters

| Analysis            | Specification   | Result |
|---------------------|---|--------|
| Functional          | A 1Kb fragment is reverse transcribed from Poly A RNA with a dilution series of BioScript™, using standard conditions. Single distinct bands were observed with agarose gel electrophoresis (ethidium stained). | Passed |
| RNase contamination | Quantitative PCR analysis with high and low RNase standards. Test sample must show less RNase activity than the limit of detection $9.7 \times 10^{-3}$ ng/μL RNase.  | Passed |

QA / QC Representative:



Andrew Galeeba-M

Date: 5<sup>th</sup> September 2022

**United Kingdom**

Tel: +44 (0)20 8830 5300  
Fax: +44 (0)20 8452 2822

**USA**

Tel: +1 901.382.8716  
Fax: +1 901.382.0027

**Germany**

Tel: +49 (0)3371 60222 00  
Fax: +49 (0)3371 60222 01

**Australia**

Tel: +61 (0)2 9209 4180  
Fax: +61 (0)2 9209 4763

**DEPC Water**

For research or further manufacturing use only

|                     |                |
|---------------------|----------------|
| Catalog No:         | BIO-65043      |
| Lot No:             | RA384-B110300  |
| Storage Conditions: | -20°C          |
| Component Lot No:   | DWT-122908A    |
| Expiry date:        | September 2024 |

**Quality Control Parameters**

| Analysis            | Specification   | Result |
|---------------------|---|--------|
| DNA contamination   | Quantitative PCR analysis with no template. Presence of <i>E. coli</i> and mouse genomic DNA checked. Test sample must amplify in line with control sample.   | Passed |
| DNase contamination | Incubation of a 1Kb double stranded DNA fragment. Incubation for 4 hours at 37°C with dilution series of DNase I. Analysed by agarose gel electrophoresis. Test sample must show less degradation than the limit of detection $2.5 \times 10^{-3}$ U DNase I. | Passed |
| RNase contamination | Quantitative PCR analysis with high and low RNase standards. Test sample must show less RNase activity than the limit of detection $9.7 \times 10^{-3}$ ng/ $\mu$ L RNase.  | Passed |

QA / QC Representative:



Andrew Galeeba-M

Date: 5<sup>th</sup> September 2022

**United Kingdom**

Tel: +44 (0)20 8830 5300  
Fax: +44 (0)20 8452 2822

**USA**

Tel: +1 901.382.8716  
Fax: +1 901.382.0027

**Germany**

Tel: +49 (0)3371 60222 00  
Fax: +49 (0)3371 60222 01

**Australia**

Tel: +61 (0)2 9209 4180  
Fax: +61 (0)2 9209 4763