MyTaq[™] One-Step RT-PCR Mix

For research or further manufacturing use only

Catalog No:	BIO-65049
Lot No:	RA387-B110430
Storage Conditions:	-20°C
Component Lot No:	MTOS-022308A
Expiry date:	September 2024

Quality Control Parameters

Analysis	Specification	Result
Functional	Fragments of sizes 1000bp and 1400bp were amplified with a dilution series of mouse RNA, using standard conditions and 45 cycles. Single distinct bands were observed with agarose gel electrophoresis (ethidium stained).	Passed
DNA contamination	Quantitative PCR analysis with no template. Presence of <i>E. coli</i> and mouse genomic DNA checked. Test sample must amplify in line with control sample.	Passed
DNase contamination	Incubation of a 1Kb double stranded DNA fragment. Incubation for 4 hours at 37°C with dilution series of DNase I. Analysed by agarose gel electrophoresis. Test sample must show less degradation than the limit of detection 2.5 x 10 ⁻³ U DNase I.	Passed
RNase contamination	Quantitative PCR analysis with high and low RNase standards. Test sample must show less RNase activity than the limit of detection $9.7x10^{-3}$ ng/µL RNase.	Passed

QA / QC Representative:

Andrew Galeeba-M

Date: 22nd August 2022

Tel: +44 (0)20 8830 5300 Fax: +44 (0)20 8452 2822 Germany Tel: +49 (0)3371 60222 00 Fax: +49 (0)3371 60222 01 Australia

Tel: +61 (0)2 9209 4180 Fax: +61 (0)2 9209 4763

info.uk@meridianlifescience.com; www.meridianlifescience.com

USA

COA No: CA_BEM-0010

Version: 07

Catalog No:BIO-65049Lot No:RA387-B110430Storage Conditions:-20°CComponent Lot No:RTS-222108AExpiry date:September 2024

Quality Control Parameters

Analysis	Specification	Result
Functional	Quantitative PCR analysis amplifying 6 genes from a dilution series of mouse RNA under standard conditions. Cq and melt profiles must be consistent for the test and reference sample with \pm 0.5 Cq variance.	Passed
DNA contamination	Quantitative PCR analysis with no template. Presence of <i>E. coli</i> and mouse genomic DNA checked. Test sample must amplify in line with control sample.	Passed
DNase contamination	Incubation of a 1Kb double stranded DNA fragment. Incubation for 4 hours at 37°C with dilution series of DNase I. Analysed by agarose gel electrophoresis. Test sample must show less degradation than the limit of detection 2.5 x 10 ⁻³ U DNase I.	Passed
RNase contamination	Quantitative PCR analysis with high and low RNase standards. Test sample must show less RNase activity than the limit of detection $9.7x10^{-3}$ ng/µL RNase.	Passed

QA / QC Representative:

Andrew Galeeba-M

Date: 22nd August 2022

Tel: +44 (0)20 8830 5300 Fax: +44 (0)20 8452 2822 <u>Germany</u> Tel: +49 (0)3371 60222 00 Fax: +49 (0)3371 60222 01

Australia

Tel: +61 (0)2 9209 4180 Fax: +61 (0)2 9209 4763



	Catalog No:	BIO-65049
RNase Inhibitor	Lot No:	RA387-B110430
	Storage Conditions:	-20°C
Suitable for Research and further Manufacturing Use	Component Lot No:	RI-122308A
	Expiry date:	September 2024

Quality Control Parameters

Analysis	Specification	Result
Inhibition	Test level of inhibition by incubating total RNA with concentration gradient of RNase A. Bands were observed with agarose gel electrophoresis (ethidium stained).	Passed

QA / QC Representative:

Andrew Galeeba-M

Date: 22nd August 2022

<u>Germany</u> Tel: +49 (0)3371 60222 00 Fax: +49 (0)3371 60222 01 <u>Australia</u>

Tel: +61 (0)2 9209 4180 Fax: +61 (0)2 9209 4763



COA No: CA_XBS-0020

Version: 07

	Catalog No:	BIO-65049
	Lot No:	RA387-B110430
DEPC Water	Storage Conditions:	-20°C
For research or further manufacturing use only	Component Lot No:	DWT-122808B
	Expiry date:	September 2024

Quality Control Parameters

Analysis	Specification	Result
DNA contamination	Quantitative PCR analysis with no template. Presence of <i>E. coli</i> and mouse genomic DNA checked. Test sample must amplify in line with control sample.	Passed
DNase contamination	Incubation of a 1Kb double stranded DNA fragment. Incubation for 4 hours at 37°C with dilution series of DNase I. Analysed by agarose gel electrophoresis. Test sample must show less degradation than the limit of detection 2.5 x 10 ⁻³ U DNase I.	Passed
RNase contamination	Quantitative PCR analysis with high and low RNase standards. Test sample must show less RNase activity than the limit of detection $9.7x10^{-3}$ ng/µL RNase.	Passed

QA / QC Representative:

Andrew Galeeba-M

Date: 22nd August 2022

<u>Australia</u>

Tel: +61 (0)2 9209 4180 Fax: +61 (0)2 9209 4763

info.uk@meridianlifescience.com; www.meridianlifescience.com